

2025-11-06



HRP-Defender (10x)

Product overview

Catalog number CR321

Description Concentrated diluent for long-term storage of HRP-conjugates.

Animal- and protein-free.

Storage 2-8 °C

pH-value at 19.0 – 21.0 °C 6.1 ± 0.2

Preservative Contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1)

Expiry date when stored

unopened

See label on the bottle

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Fields of application

Dilution of Horseradish Peroxidase (HRP)-conjugated antibodies, proteins and biomolecules, as well as free HRP for long-term storage in solution at non-freezing temperatures, ensuring stability and functionality for subsequent use in immunoassays and related analytical applications.

Instructions for use

HRP-Defender (10x) dilution

To prepare a working solution, HRP-Defender (10x) can be diluted up to 1:10 into another buffer - e.g. Sample Buffer (CR105 / CR305) - before use

(1 part of HRP-Defender + 9 parts of buffer).

Please shake HRP-Defender (10x) thoroughly prior to dilution.

HRP-Defender (10x) must **not** be diluted into **water**.

For unusually unstable HRP-conjugates, HRP-Defender (10x) may be used at a higher concentration - such as a 1:5 dilution - or even undiluted.

Dilution of the conjugate

HRP-conjugates can be diluted and stored directly in the working solution of HRP-Defender. Typical conjugate concentrations in the final solution range from 40 to 500 ng/ml. Once diluted in a suitable buffer, HRP-Defender can

directly be used as assay buffer in immunoassays.

Storage of the conjugate

HRP-Defender (10x) is designed for the storage of conjugates in solution at non-freezing temperatures (i.e. above 0 °C).

Stability data of one peroxidase conjugate cannot be directly transferred to other conjugates. Each conjugate must be tested for its individual shelf-life in the respective HRP-Defender (10x) dilution. If HRP-Defender (10x) is

Page 1 of 2 V-11-2025



CANDOR Bioscience GmbH,



Product specification STABILIZER

2025-11-10

used in immunodiagnostic kits, the shelf life must be validated according to the applicable regulatory requirements for diagnostics.

HRP-Defender (10x) contains components that may interfere with commonly used conjugation methods, such as techniques targeting primary amines. Suitability for any given conjugation method must be tested in advance. We recommend diluting biomolecules in HRP-Defender dilutions only **after** conjugation.

Assay interference

If background or false-positive signals occur - e.g. due to cross-reactivities, matrix effects, or HAMA - we recommend diluting HRP-Defender (10x) into LowCross-Buffer® (CR100) or LowCross-Buffer® protein-free (CR300). Suitability of HRP-Defender for a specific assay must be tested by the user.

Stability

HRP-Defender (10x) contains the preservative CMIT/MIT but is **not** sterile. To guarantee optimal stability, the respective dilution buffer should contain appropriate preservative levels.

Please note that high protein concentrations and/or microbial contamination may reduce the effectiveness of the preservative. If antibodies/conjugates are added for storage in a non-sterile manner and microbial contamination is suspected, higher-than-usual preservative concentrations may be required.

LowCross is a registered trade mark of CANDOR Bioscience.



CANDOR Bioscience GmbH,
part of Medix Biochemica
Simoniusstrasse 39
88239 Wangen, Germany

Page 2 of 2

medix@medixbiochemica.com www.medixbiochemica.com

V-11-2025