

HiDi® 2x PCR Master Mix

Introduction

HiDi® stands for High Discrimination of nucleotides mismatching at the 3'-terminus of PCR primers. The HiDi® 2X PCR Master Mix is a SYBR Green-compatible ready to use mix that contains a highly selective recombinant DNA polymerase suitable for allele-specific PCR (AS-PCR)¹, primer extension or methylation-specific PCR.

HiDi® is similarly effective for checking the quality of genome editing techniques such as CRISPR/Cas or TALEN-based approaches^{2,3}. The enzyme uses a hot-start aptamer-based

technology while the optimized master mix chemistry promotes high sensitivity and robust polymerase processivity. Temperatures above 50-55°C melt the aptamer's secondary structure and set free the DNA polymerase.

HiDi® efficiently discriminates transversion mutations as well as transition mutations in a more delayed fashion. For probe-based assays, HiDi® Taq 2X PCR Master Mix is also available (Cat no. #4200S or #4200M).

HiDi® efficiently detects mutations by discriminating wild type mismatching primers, which is resulting in diminished or completely absent amplification.



Methods

HiDi® 2x PCR Master Mix and three competitors were used to amplify a 100 bp-long fragment of actin gene (Genbank NM_001101.5). A 10-fold dilution series of human genomic DNA was amplified in triplicate using SYBR Green chemistry over 40 cycles as per

Table 1. Matching versus 3'-mismatching reverse primers were designed according to Figure 1. Reaction setup and final primer concentrations were applied according to the manufacturer recommendations.

Step	Temperature	Time	Cycles
1	95°C	2 min	1
2	95°C	15 sec	40
3	60°C	30 sec	40
4	72°C	30 sec	40

Table 1: Thermal cycling conditions. Step 1 of initial denaturation was extended to 10 min for competitor C.

NM_001101.5:85-1212	3' -GTCTTTGCGGATGTCCACGTCACACTTCATGATGGAGT-5'
WT primer	5' -GGATGTCCACGTCACACTTC-3'
A-mismatch primer	5' -GGATGTCCACGTCACACTTA-3'
T-mismatch primer	5' -GGATGTCCACGTCACACTTT-3'
Forward primer	5' -CACTCTCCAGCCTTCCTTC-3'

Figure 1: Primer design. The mismatching reverse primers exhibit either a transversion (A-mismatch primer) or transition (T-mismatch primer) mutation vs. the matching WT primer. All primers have a Tm of ~63°C. WT stands for Wild Type.

HiDi® efficiently detects the wild type only by strong and early amplifications. T-allele is absent as the amplification reaction with the T-allelic primer is significantly delayed and hence demonstrates good discrimination.

Results

Firstly, HiDi® 2x PCR Master Mix outperformed and was the most efficient for discriminating against a transversion mutation (A-mismatch primer) in comparison with the WT primer as shown in panels A and E of Figure 2.

The product of the T-mismatch primer is significantly delayed and hence, demonstrates good discrimination. Furthermore, competitors' products

yielded lower fluorescence signals and resulted in poor discrimination of either the A-mismatch or T-mismatch primer as shown in panels B, C, D and E of Figure 2.

Secondly, the efficiency of HiDi® 2x PCR Master Mix after amplification with the matching WT primer stands at 97% and proves to be superior and over the three other products tested for three major suppliers.

Conclusions

- HiDi® 2x PCR Master Mix displays superior capacity for discriminating against a single mutated nucleotide. It is a product of choice for mutation detection using AS-PCR and for sensitive genome editing techniques.
- HiDi® 2x PCR Master Mix demonstrates superior performance with SYBR Green chemistry compared to three mixes from leading suppliers.
- HiDi® 2x PCR Master Mix is an equally efficient mix for highly sensitive PCR amplification.

HiDi® 2x PCR Master Mix is the master mix of choice for single nucleotide detection by allele-specific PCR.




Figure 2: Amplification curves and Cq values plots.
A) HiDi® 2x PCR Master Mix discriminates primers with a 3'-end mismatch.
B) Competitor A does not discriminate either of the mismatches.
C) Competitor B performances result in poor discrimination.
D) Competitor C performances result in poorer discrimination and signals.

E) The Cq values for HiDi® 2x PCR Master Mix and competitors across DNA copies demonstrates the high discriminative power of HiDi® in presence of a transversion or transition mutation.
 in panels A, B, C and D: All three amplification curves presented for WT, A-mismatch and G-mismatch primers correspond to the 50000, 5000 and 500 DNA copy numbers in the samples.



References:

1. Miotto, O. et al. Emergence of artemisinin-resistant *Plasmodium falciparum* with kelch13 C580Y mutations on the island of New Guinea. *PLoS Pathog.* 16, 1–21 (2020).
2. Morisaka, H. et al. CRISPR-Cas3 induces broad and unidirectional genome editing in human cells. *Nat. Commun.* 10, (2019).
3. Sakurai, T. et al. Bindel-PCR: a novel and convenient method for identifying CRISPR/Cas9-induced biallelic mutants through modified PCR using *Thermus aquaticus* DNA polymerase. *Sci. Rep.* 9, 1–14 (2019).

Medix Biochemica

More information on our products
[medixbiochemica.com](https://www.medixbiochemica.com)

Klovinpellontie 3, FI-02180 Espoo, Finland

Orders and product information
mdx@medixbiochemica.com



myPOLS Biotec GmbH
Byk-Gulden-Str. 2 // 10, 78467 Konstanz,
Germany

Copyright © 6/2022 Medix Biochemica. All rights reserved.

Medix Biochemica reserves the right to make changes and improvements to any of the products described in this document without prior notice.