## External Controls for Actim Pancreatitis test, Instructions for use

### Intended use

Actim Pancreatitis Controls are intended to be used with the Actim Pancreatitis test (32701ETAC) for external quality control. The controls may also be used to demonstrate negative results and weak and strong positive results.

### Kit components

Actim Pancreatitis Control kit (32700ETAC) consists of one of each of the following components:

- Actim Pancreatitis Negative Control (32780ETAC)
- Actim Pancreatitis Positive Control, Low (32781ETAC)
- Actim Pancreatitis Positive Control, High (32782ETAC)
- Actim Reconstitution Solution (30889ETAC)
- Instructions for use

A pipette capable to measure 500 µl of solution is needed for accurate dispensing of Reconstitution Solution (not provided with the kit).

### Storage

Store the components at +2...+8 °C. Stored unopened, each component can be used until the expiry date marked on the component. After reconstitution, the controls should be used within 24 hours if they are kept at room temperature (+18...+25 °C) or within 6 days if they are kept at +2...+8 °C.

### Control material

The positive controls consist of purified human trypsinogen-2 in buffered solution containing protein as stabilizer, protease inhibitors and <0.1% sodium azide (NaN<sub>3</sub>). The negative control consists of the same matrix without added antigen.

### Limitations of the test

- The controls are to be used with the Actim Pancreatitis test (manufactured by Medix Biochemica) only.
- The controls are intended for professional *in vitro* use only.

### Test procedure and interpretation of the results

#### Reconstitution:

1. Unscrew the cap and open the vials by carefully removing the grey stopper.
2. Reconstitute the controls by adding 500 µl of Reconstitution Solution to the vials. When pipetting, avoid the risk of carryover by changing the pipette or its tip after each addition.
3. Close the stopper and allow the vial to reconstitute for at least 15 minutes.
4. Mix the solution.

#### Testing:

1. The reconstituted controls are tested by dipping the dipsticks directly into the control vials. Each dipstick requires approximately 150 µl of sample. Before dipping, ensure that the volume of the controls in the vial is sufficient.
2. The testing is performed by carefully following the instructions for use of the Actim Pancreatitis test.

The test results are interpreted as instructed in the paragraph "Test procedure and interpretation of results" of the Actim Pancreatitis test:

- The Negative Control gives a negative test result.
- The Low Positive Control gives a weak positive result.
- The High Positive Control gives a strong positive result.

### Notes

- Before use, ensure that the package and vials are intact. Do not use a control if the vial is damaged.
- Do not use the controls before they have been completely reconstituted.
- When testing, the controls should be at room temperature to ensure optimal test results with Actim Pancreatitis test.
- Carefully follow the instructions for use of the Actim Pancreatitis test and Actim Pancreatitis Controls to ensure the correct results.
- If the results of the controls are not appropriate, repeat the procedure with another dipstick.
- The frequency of quality control testing is to be determined by each laboratory independently, and will depend on factors such as patient sample frequency and/or local regulations.
- The trypsinogen-2 antigen in the Positive Controls has been shown to be negative for HBV-DNA, HIV-RNA and HCV-RNA. Nevertheless, such tests are unable to prove the complete absence of viruses, and the controls should therefore be handled taking appropriate precautions.
- Like all biological materials, the controls should be treated as potentially hazardous, and disposed of in accordance with local authority guidelines.

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