## External controls for Actim® Calprotectin test, Instructions for use

### INTENDED USE

Actim® Calprotectin Controls are intended to be used with the Actim® Calprotectin test (35032ETAC) for external quality control. The controls may also be used to demonstrate negative results and weak and strong positive results.

### KIT COMPONENTS

Actim Calprotectin Controls kit (35000ETAC) consists of one of each of the following components:

- Actim Calprotectin Negative Control (35080ETAC)
- Actim Calprotectin Positive Control, Low (35081ETAC)
- Actim Calprotectin Positive Control, High (35082ETAC)
- Actim Reconstitution Solution (30889ETAC)
- Instructions for use

A pipette capable to measure 500 µl of solution is needed for accurate dispensing of Reconstitution Solution (not provided with the kit). For dilution of the controls, a pipette capable to measure 100 µl is also needed.

### STORAGE

Store the components at +2...+8 °C. Stored unopened, each component can be used until the expiry date marked on the component. After reconstitution, the controls should be used within 7 days. Alternatively, the reconstituted controls can also be stored frozen at -20 °C for 2 months. Do not freeze after thawing.

### CONTROL MATERIAL

The Positive Control consist of purified human calprotectin in buffered solution containing protein as stabilizer and <0.1 % sodium azide (NaN<sub>3</sub>). The Negative Control consists of the same matrix without added antigen.

### LIMITATIONS OF THE TEST

- The controls are to be used with the Actim Calprotectin test (manufactured by Medix Biochemica) only.
- The controls are intended for professional *in vitro* use only.

### TEST PROCEDURE AND INTERPRETATION OF THE RESULTS

#### Reconstitution:

1. Unscrew the cap and open the vials by carefully removing the grey stopper.
2. Reconstitute the controls by adding 500 µl of Reconstitution Solution to the vials. When pipetting, avoid the risk of carryover by changing the pipette or its tip after each addition.
3. Close the stopper and allow the vial to reconstitute for at least 15 minutes.
4. Mix the solutions.

#### Testing:

1. Dilute each control by adding 100 µl of the control to the Actim Calprotectin Specimen Dilution Buffer (35071ETAC). When pipetting, avoid the risk of carryover by changing the pipette or its tip after each addition.
2. Testing is performed by carefully following the instructions for use of the Actim Calprotectin test.

The test results are interpreted as instructed in the paragraph "Test procedure and interpretation of results" of the Actim Calprotectin test:

- The negative control gives a negative test result.
- The low positive control gives positive result with test line 1 and negative with test line 2
- The high positive control gives positive result with both test line 1 and 2.

### NOTES

- Before use, ensure that the package and vials are intact. Do not use a control if the vial is damaged.
- Do not use the controls before they have been completely reconstituted.
- When testing, the controls should be at room temperature to ensure optimal test results with Actim Calprotectin test.
- Carefully follow the instructions for use of the Actim Calprotectin test and Actim Calprotectin Controls to ensure correct results.
- If the results of the controls are not appropriate, repeat the procedure with another dipstick.
- The frequency of quality control testing is to be determined by each laboratory independently, and will depend on factors such as patient sample frequency and/or local regulations.
- Like all biological materials, the controls should be treated as potentially hazardous, and disposed of in accordance with local authority guidelines.

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