

---

## Product Datasheet

### PlateBlock™ CR112

<b>Product Name</b>	PlateBlock™
<b>Catalog Number</b>	CR112
<b>Description</b>	Animal-free and protein-free blocker for saturating free binding sites on plastic surfaces, optimized for reducing background in serological assays
<b>Tested Applications</b>	ELISA
<b>Alternative Names</b>	surface blocker
<b>Brand</b>	CANDOR
<b>Storage</b>	2 - 8 °C or -15 to -30 °C (tolerates repeated freezing and thawing cycles)
<b>Note</b>	PlateBlock™ is ready-to-use. Shake thoroughly before use
<b>Preservatives</b>	Contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1)