

Product Datasheet

Washing Buffer PBS 10x CR140

Product Name	Washing Buffer PBS 10x
Catalog Number	CR140
Description	Wash Buffer for immunoassays
Tested Applications	WB, IHC, Multiplex Assay, Immuno-PCR, ELISA
Brand	CANDOR
Storage	2 – 8 °C or -15 to -30 °C (tolerates repeated freezing and thawing cycles)
Note	Salt crystals may precipitate due to storage at 2 - 8 °C or freezing. Therefore, to prepare the working solution, the Washing Buffer must be warmed to room temperature beforehand, which will dissolve any salts that may have precipitated. The buffer should be thoroughly mixed again by shaking immediately before use. The working solution is prepared from the stock solution by diluting 1:10 with distilled or demineralized water and should be used on the same day. Washing Buffer can be used in automated washers depending on the salt tolerance of the washer. It is essential to observe the specifications of the respective manufacturer! Suitability of Washing Buffer for a specific assay has to be tested by the user

Preservatives Contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1)