

Oy Medix Biochemica Ab

Headquarters Klovinpellontie 3 FI-02180 Espoo

Finland

Email: sales@medixbiochemica.com

Medix Biochemica USA Inc.

10850 Metro Court Maryland Heights, MO 63043 United States of America

Email: sales@medixbiochemica.com

Product Datasheet

HiFi DNA Polymerase Hot-Start #3000

Product Name

HiFi DNA Polymerase Hot-Start

Catalog Number

#3000

Description

HiFi DNA Polymerase Hot-Start is an aptamer-based hot-start, high fidelity hyperthermophilic recombinant DNA polymerase from the archaeon Pyrococcus furiosus. HiFi DNA Polymerase Hot-Start exhibits 5' to 3' polymerase activity and 3' to 5' exonuclease activity. The optimized buffer chemistry facilitates high sensitivity, yield, specificity, robust and rapid polymerase. processivity. The enzyme is ideal for long, complex, difficult DNA templates and is resistant to PCR inhibitors. HiFi DNA Polymerase Hot-Start has a lower error rate than standard Taq DNA polymerase (100x) and therefore, is suitable for PCR applications

where higher accuracy is needed.

Tested Applications

End-Point, Real-Time

Brand Medix Biochemica

Storage -20°C

Content

100 U: 2 U/µl, 1 x 0.05 ml HiFi DNA Polymerase Hot-Start; 1 x 1.7

ml 5x HiFi Buffer; 1 x 1.7 ml 10x HiFiOpt Enhancer

500 U: 2 U/µl, 1 x 0.25 ml HiFi DNA Polymerase Hot-Start; 3 x 1.7

ml 5x HiFi Buffer; 2 x 1.7 ml 10x HiFiOpt Enhancer