



Product specifications

Name	Anti-h cTnI 9707 SPTN-5	
Specificity	Antibody recognizes human cardiac troponin I (cTnI)	
Description	Monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components	
Product code	100180	
Product buffer solution	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN $_{3}$ as a preservative	
Shelf life and storage	18 months from manufacturing at 2–8 °C	
Subclass	IgG ₁	
Analyte description	Troponins form a complex of three regulatory proteins (troponin C, I and T) that are integral to muscle contraction in skeletal and cardiac muscles. Serum cardiac troponin tests can be used to help diagnose several different heart disorders, especially myocardial infarction.	

Parameters tested on each lot

Product appearance	Liquid, may turn slightly opaque during storage	
Product concentration	5.0 mg/ml (+/- 10 %)	
Immunoreactivity	80–120 % compared to the reference sample in an FIA test	
IEF Profile	5.9–7.4	
Purity	≥ 95 %	

Kinetic parameters

Association rate constant	1.9 x 10 ⁶ 1/Ms
Dissociation rate constant	5.8 x 10 ⁻⁶ 1/s
Affinity constant	KA = 3.3 x 10 ¹¹ 1/M; KD = 3.0 x 10 ⁻¹² M (= 0.003 nM)
Determination method	SPR analysis (ProteOn XPR36)
Determination antigen	Human cardiac troponin I-T-C complex, HyTest (Cat 8T62, Lot 13/07)

Oy Medix Biochemica Ab Klovinpellontie 3 FI-02180 Espoo, Finland

medix@medixbiochemica.com www.medixbiochemica.com VAT reg.no. FI14631532

Legal disclaimer Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. Medix/MAB and Medix/Antigens are trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.





Cross-reactivities

Skeletal troponin I (sTnI) 14 %

Does not recognize cardiac troponin T (cTnT), cardiac troponin C (cTnC), or skeletal troponin T (sTnT)

Epitope

In the C-terminal part of the cTnI molecule (residues 190-195)

Pair recommendations

DETECTION 9701/ 9707/ 9703 9705 9708 9709 RC9750 RC9701 RC9707 9701/ + + + + + + _ RC9701 9703 _ _ + _ + + _ 9705 CAPTURE + + _ + --+ 9707/ RC9707 9708 -_ ---_ _ 9709 _ + _ _ _ + _ RC9750 _ + + + _ + +

RC9750 is recommended as a very sensitive capture.

For increased sensitivity, 2+1 or 2+2 antibody approach is recommended, for example capture: 9701/RC9701 & 9705, detection 9703 & 9707/RC9707.

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Platforms tested	FIA			
Antigens tested		Recombinant cTnI antigen, Medix Biochemica 610102, native cTnI antigen Lee Biosolutions 550-11, and native cardiac Troponin Complex (cTn ITC) antigen, Lee Biosolutions 550-08.		
Product stability	TEMPERATURE, TIME	RESULT		
	-70 °C, 21 days	ОК		
	-20 °C, 21 days	ОК		
	+4 °C, 21 days	ОК		
	+35 °C, 21 days	ОК		
	+45 °C, 7 days	ОК		
	Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.			
Miscellaneous	-			
References	Savukoski, T. et al. (2012) Troponin-specific autoantibody interference in different cardiac troponin I assay configurations. Clin. Chem. 58:1040-1048.			

Oy Medix Biochemica Ab Klovinpellontie 3 FI-02180 Espoo, Finland medix@medixbiochemica.com www.medixbiochemica.com VAT reg.no. FI14631532 Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. Medix/MAB and MedixAntigens are trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.