



Product specifications

Name	Anti-h cTnl RC9701 SPTN-5			
Specificity	Antibody recognizes human cardiac troponin I (cTnI)			
Description	Recombinant chimeric antibody containing variable region from mouse IgG and constant region from human IgG ₁ , cultured <i>in vitro</i> under conditions free from animal-derived components.			
Product code	140000			
Product buffer solution	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN $_3$ as a preservative			
Shelf life and storage	Unspecified, storage at 2–8 °C			
Subclass	human IgG ₁			
Analyte description	Troponins form a complex of three regulatory proteins (troponin C, I and T) that are integral to muscle contraction in skeletal and cardiac muscles. Serum cardiac troponin tests can be used to help diagnose several different heart disorders, especially myocardial infarction.			

Parameters tested on each lot

Product appearance	Liquid, may turn slightly opaque during storage
Product concentration	5.0 mg/ml (+/- 10 %) (A280 nm, 1 mg/ml, 1 cm = 1.3)
Immunoreactivity	80–120 % compared to the reference sample in an FIA test
IEF Profile	7.5–8.8
Purity	≥ 95 %

Kinetic parameters

Association rate constant	Not Determined (N/D)
Dissociation rate constant	N/D
Affinity constant	N/D
Determination method	-
Determination antigen	-

medix@medixbiochemica.com www.medixbiochemica.com VAT reg.no. FI14631532

Legal disclaimer Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. Medix/MAB and Medix/Antigens are trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.





Cross-reactivities

Epitope

Does not recognize skeletal troponin I (sTnI).

In the stable part of the cTnI molecule (residues 85-95)

Pair recommendations

		DETECTION						
		9701/ RC9701	9703	9705	9707/ RC9707	9708	9709	RC9750
CAPTURE	9701/ RC9701	-	+	+	+	+	+	+
	9703	-	-	+	-	+	+	-
	9705	+	+	-	+	-	-	+
	9707/ RC9707	-	-	-	-	-	-	-
	9708	-	-	-	-	-	-	-
	9709	-	+	-	-	-	-	+
	RC9750	+	-	+	+	+	+	-

RC9750 is recommended as a very sensitive capture.

For increased sensitivity, 2+1 or 2+2 antibody approach is recommended, for example capture: 9701/RC9701 & 9705, detection 9703 & 9707/RC9707.

Please note that pair recommendations are based on results obtained by our laboratory using cardiac troponin complex (I-C) as an antigen. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Platforms tested

Antigens tested

Product stability

FIA

Recombinant cTnI antigen, Medix Biochemica 610102, native cTnI antigen Lee Biosolutions 550-11, and native cardiac Troponin Complex (cTn ITC) antigen, Lee Biosolutions 550-08.

TEMPERATURE, TIME	RESULT
-70 °C, 21 days	ОК
-20 °C, 21 days	ОК
+4 °C, 21 days	ОК
+35 °C, 21 days	ОК
+45 °C, 7 days	ОК

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.

Miscellaneous

References

Oy Medix Biochemica Ab Klovinpellontie 3 FI-02180 Espoo, Finland medix@medixbiochemica.com www.medixbiochemica.com VAT reg.no. FI14631532 Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. Medix/MAB and MedixAntigens are trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.