# Medix Biochemica

### Anti-h IgE 8510 SPTN-5

Product overview		
Catalog number	100429	
Specificity	Antibody recognizes human immunoglobulin E	
Description	Monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components.	
Product buffer solution	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN $_3$ as a preservative	
Shelf life and storage	24 months from manufacturing at 2–8 °C	
Subclass	lgG <sub>1</sub>	
Analyte description	Immunoglobulin E (IgE) is a class of antibody that has only been found in mammals. It plays an important role in allergy and defense against parasites. Although IgE is typically the least abundant isotype - blood serum IgE levels in a normal individual are only 0.05 % of the IgG concentration - it is capable of triggering the most powerful immune reactions.	
Parameters tested on ea	ach lot	
Product appearance	Liquid, may turn slightly opaque during storage	
Product concentration	5.0 mg/ml (+/- 10 %)	
Immunoreactivity	80–120 % compared to the reference sample in an FIA test	
IEF Profile	5.8–7.0	
Purity	≥ 95 %	
Kinetic parameters		
Association rate constant	Not Determined (N/D)	
Dissociation rate constant	N/D	
Affinity constant	K <sub>A</sub> = 8 x 10 <sup>9</sup> 1/M	
Determination method	SPR analysis (Biacore)	
Determination antigen	IgE, Biodesign (Cat A101644, Lot 912737)	



**Oy Medix Biochemica Ab** Klovinpellontie 3, FI-02180 Espoo, Finland medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

#### Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.

## Medix Biochemica

Cross-reactivities	Human IgA	<0.05 %
	Human IgG	<0.05 %
	Human IgM	<0.05 %

N/D

**Epitope** 

### **Pair recommendations**

		DETECTION		
		8510	8516	8518
CAPTURE	8510	-	+	+
	8516	-	-	-
	8518	-	-	-

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Platforms	tested	FIA	

Antigens tested N/D

Prod

duct stability	TEMPERATURE, TIME	RESULT
-	-70 °C, 21 days	N/D
	-20 °C, 21 days	OK
	+4 °C, 21 days	OK
	+30 °C, 7 days	OK
	+35 °C, 21 days	OK
	+45 °C, 7 days	OK

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.

**Miscellaneous** 

References



**Oy Medix Biochemica Ab** Klovinpellontie 3, FI-02180 Espoo, Finland

medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

#### Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.