# Medix Biochemica

### Anti-h PG II 8102 SPRN-1

Product overview					
Catalog number	100098				
Specificity	Antibody recognizes human pepsinogen II				
Description	Monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components.				
Product buffer solution	37 mM citrate, 125 mM phosphate, pH 6.0, 0.9 % NaCl, 0.095 % NaN $_{\rm 3}$ as a preservative				
Shelf life and storage	24 months from manufacturing at 2–8 °C				
Subclass	lgG1				
Analyte description	Pepsinogen is the pro-form of pepsin and is produced in the stomach by chief cells. The major part of pepsinogen is secreted into the gastric lumen but a small amount can be found in the blood. Alterations in the serum pepsinogen concentrations has been found with <i>Helicobacter pylori</i> ( <i>H. pylori</i> ) infections, peptic ulcer disease, gastritis, and gastric cancer. More precise analysis may be achieved by measuring the pepsinogen I/II ratio.				
Parameters tested on each lot					
Product appearance	Liquid, may turn slightly opaque during storage				
Product concentration	1.0 mg/ml (+/- 10 %)				
Immunoreactivity	80–120 % compared to the reference sample in an FIA test				
IEF Profile	6.2–8.0				
Purity	≥ 95 %				
Kinetic parameters					
Association rate constant	Not Determined (N/D)				
Dissociation rate constant	N/D				
Affinity constant	1 x 10 <sup>9</sup> 1/M				
Determination method	Radioimmunoassay (RIA)				
Determination antigen	PG II, Medix Biochemica in-house (A87/94)				



**Oy Medix Biochemica Ab** Klovinpellontie 3, Fl-02180 Espoo, Finland medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

### Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.

# Medix Biochemica

Cross-reactivities	Human Pepsinogen I Porcine Pepsinogen Gastrin I Somatostatin Peptide YY Bombesin			0.9 % (Medix Biochemica in-house, A125/92) < 0.22 % (Sigma, Cat P-4656, Lot 82H8065) < 0.22 % (Sigma, Cat G-9020, Lot 110H48131) < 0.22 % (Sigma, Cat S-9129, Lot 81H00471) < 0.22 % (Sigma, Cat P-1306, Lot 90H48992) < 0.22 % (Sigma, Cat B-4272, Lot 82H58501)				
Epitope	N/D							
Pair recommendations	-		DETECTION					
			8101	8102	8103	_		
	CAPTURE	8101	-	+	+	-		
		8102	-	-	+	_		
	Q A	8103	+	+	-			
	<ul> <li>Following pairs are especially recommended for the below mentioned assays:</li> <li>FIA: 8101 (capture) – 8102 (detection), 8101 – 8103 and 8103 – 8102 LF: 8103 (membrane) – 8101 (particles) and 8103 – 8102 IT: 8101 – 8103</li> <li>Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.</li> </ul>							
Platforms tested	FIA, LF, IT							
Antigens tested	Recombinant Pepsinogen II antigen, Medix Biochemica 610001 and native Pepsinogen II antigen Lee Biosolutions 440-52.							
Product stability	TEMPERATURE, TIME -70 °C, 21 days -20 °C, 21 days +4 °C, 21 days +25 °C, 21 days +35 °C, 21 days +35 °C, 21 days +45 °C, 7 days Stability testing is performed in the temperatures affect the antigen bind antibody. Please note that the shelf real time stability testing at 2–8 °C i				ng, charge fe given on	or composition of the the first page is based on		
Miscellaneous	-							
References	-							



**Oy Medix Biochemica Ab** Klovinpellontie 3, FI-02180 Espoo, Finland

medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

Legal disclaimer Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.