Medix Biochemica

Anti-h Thyroglobulin 2803 SPTN-5

Product overview					
Catalog number	100332				
Specificity	Antibody recognizes human thyroglobulin				
Description	Monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components.				
Product buffer solution	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN $_3$ as a preservative				
Shelf life and storage	36 months from manufacturing at 2–8 °C				
Subclass	lgG _{2b}				
Analyte description	Thyroglobulin is a 660 kDa dimeric protein produced by and used entirely within the thyroid gland to produce the thyroid hormones thyroxine (T4) and triiodothyronine (T3). Thyroglobulin levels in the blood can be used as a tumor marker for certain kinds of thyroid cancer. Thyroglobulin levels in the blood can also be elevated in cases of Graves' disease.				
Parameters tested on each lot					
Product appearance	Liquid, may turn slightly opaque during storage				
Product concentration	5.0 mg/ml (+/- 10 %)				
Immunoreactivity	80–120 % compared to the reference sample in an FIA test				
IEF Profile	6.5–8.0				
Purity	≥ 95 %				
Kinetic parameters					
Association rate constant	1.7 x 10⁵ 1/Ms				
Dissociation rate constant	3.7 x 10 ⁻⁵ 1/s				
Affinity constant	$K_A = 4.5 \times 10^9 \text{ 1/M}; K_D = 2.2 \times 10^{-10} \text{ M} (=0.22 \text{ nM})$				
Determination method	SPR analysis (ProteOn XPR36)				
Determination antigen	Thyroglobulin, Fitzgerald (Cat 30-AT01, Lot A11102002)				



Oy Medix Biochemica Ab Klovinpellontie 3, FI-02180 Espoo, Finland medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.

Medix Biochemica

Cross-reactivities

Does not recognize thyroid peroxidase, others not tested.

Epitope

Pair recommendations

		DETECTION				
		2802	2803	2804	2805	
CAPTURE	2802	-	-	+	+	
	2803	-	-	+	+	
	2804	-	-	-	-	
0	2805	+	+	+	-	

Following pairs are especially recommended for the below mentioned assays: FIA: 2802 (capture) – 2805 (detection)

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Platforms tested FIA

Antigens tested Native Thyroglobulin antigen, Lee Biosolutions 528-11

Not Determined (N/D)

Product stability

t stability	TEMPERATURE, TIME	RESULT	
	-70 °C, 21 days	OK	
	-20 °C, 21 days	OK	
	+4 °C, 21 days	OK	
	+25 °C, 21 days	OK	
	+35 °C, 21 days	OK	
	+45 °C, 7 days	OK	

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 $^{\circ}$ C in the product buffer.

Miscellaneous

T3 or T4 has no effect on the binding of this antibody.

References



Oy Medix Biochemica Ab Klovinpellontie 3, Fl-02180 Espoo, Finland

medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.