Medix Biochemica

Anti-h cTnT 3710 SPTN-5

Product overview				
Catalog number	100698			
Specificity	Antibody recognizes human cardiac troponin T (cTnT)			
Description	Monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components.			
Product buffer solution	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN $_3$ as a preservative			
Shelf life and storage	Unspecified, storage at 2–8 °C			
Subclass	lgG ₁			
Analyte description	Troponins form a complex of three regulatory proteins (troponin C, I and T) that are integral to muscle contraction in skeletal and cardiac muscles. Serum cardiac troponin tests can be used to help diagnose several different heart disorders, especially myocardial infarction.			
Parameters tested on ea	ch lot			
Product appearance	Liquid, may turn slightly opaque during storage			
Product concentration	5.0 mg/ml (+/- 10 %)			
Immunoreactivity	80–120 % compared to the reference sample in an FIA test			
IEF Profile	6.5–7.0			
Purity	≥ 95 %			
Kinetic parameters				
Association rate constant	Not Detemined (N/D)			
Dissociation rate constant	N/D			
Affinity constant	N/D			
Determination method	-			
Determination antigen	-			



Oy Medix Biochemica Ab Klovinpellontie 3, Fl-02180 Espoo, Finland medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.

Medix Biochemica

Cross-reactivities

Does not recognize skeletal troponin T (sTnT) or cardiac troponin I (cTnI).

Epitope

Amino acids 139-152 (RIRNEREKERQNRL) of cardiac troponin T isoform 6.

Pair recommendations

		DETECTION							
		3701	3703	3708	3710	3711	3712		
CAPTURE	3701	-	-	-	-	-	-		
	3703	-	-	-	+	+	+		
	3708	+	+	-	+	+	+		
	3710	-	+	+	-	-	-		
	3711	-	+	-	-	-	-		
	3712	-	+	+	-	-	-		

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Platforms tested

FIA

+4 °C, 21 days

+45 °C, 7 days

+35 °C, 21 days

Antigens tested Recombinant cardiac troponin T antigen, Medix Biochemica 610101 and native cardiac troponin T antigen, Lee Biosolutions 550-21.

Product stability TEMPERATURE, TIME -70 °C, 21 days -20 °C, 21 days RESULT OK OK OK OK OK

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.

Miscellaneous

References



Oy Medix Biochemica Ab Klovinpellontie 3, Fl-02180 Espoo, Finland

medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.