

Actim<sup>®</sup>1ngeni

# Ingenious – here's why

- **Automation** saves time and ensures consistency
- **Compact size** requires little space
- **Simple interface** is logical and easy to navigate
- **Traceable data** provides clear, comprehensive documentation
- **Customizable** to your specific professional needs

## At the heart of the Actim 1ngeni System

The Actim 1ngeni instrument is the heart of the Actim 1ngeni system. It automatically analyses, displays and saves results of Actim PROM and Partus 1ngeni tests.

- Actim PROM 1ngeni is a reliable test for detecting premature rupture of fetal membranes that can be used for all patients.
- Actim Partus 1ngeni test can quickly determine the risk of preterm delivery.

## Key benefits

- Allows focus to be on patient, not testing procedures.
- Test results are clear, comprehensive and easy to trace – for improved patient security.
- Test results are interpreted digitally – so it takes less time, and no manual interpretation of results is needed.

Actim<sup>®</sup>1ngeni



## Ordering information

Actim 1ngeni instrument, ref. 19100AC

## Test kits to be ordered separately:

Actim PROM 1ngeni 10 tests kit, ref. 30831RETAC

Actim Partus 1ngeni 10 tests kit, ref. 31931RETAC

## Add-ons

USB-connected printer, barcode reader/scanner, external battery

## Medix Biochemica

Tel. +358 9 547 680  
Fax +358 9 505 3441  
medix@medixbiochemica.com

© 03/2015 Medix Biochemica.  
All rights reserved. Actim is a registered trademark of Medix Biochemica.

 [www.medixbiochemica.com](http://www.medixbiochemica.com)

 Medix  
Biochemica

Actim<sup>®</sup>1ngeni

# The Actim 1ngeni System

**DIGITAL CONFIDENCE IN PREGNANCY MONITORING**

The Actim 1ngeni system takes testing to a whole new level. Testing is easy, fast and automated – delivering a high degree of confidence and clarity to pregnancy monitoring.

Very accurate and user-friendly, Actim 1ngeni makes testing and documentation worry-free and reliable – both in the laboratory and at the point of care.

 Medix  
Biochemica

# Actim 1ngeni is... simply ingenious!



## Automated and traceable

Automated test results are available within five minutes. No supervision, interpretation or transcription of results is required, so doctors and midwives can devote more time to their patients. The instrument automatically makes sure that the performed test is valid – and to ensure confidence in the results, recorded data is traceable.

## Consistent and reliable results

The instrument automatically interprets the results so there are no variations in results between users.

Several QC options provide full control over the system; the test and instrument validity can be monitored via liquid controls, self-tests and instrument checks. Together, these features ensure consistent reliability in results.


## Customizable to your needs

The instrument functions and QC procedure are customizable to the requirements of your health care facility. E.g. the device can be set to remind users at prescribed intervals to do the QC. Bottom line is a worry-free work process.



**SAVE YOUR DIGITIZED PROM AND PARTUS TEST RESULTS WITHIN JUST FIVE MINUTES**

## Six easy steps:



## Results shown clearly:

- PROM positive > Membranes are ruptured
- PROM negative > No membrane rupture
- Partus positive > Elevated risk of preterm delivery
- Partus negative > Preterm delivery highly unlikely