

PRODUCT SPECIFICATIONS

Name	Anti-h NSE 9602 SPTN-5
Specificity	Antibody recognizes human neuron-specific enolase, γ -isoform
Description	Monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components
Product code	100408
Product buffer solution	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
Shelf life and storage	24 months from manufacturing at 2-8 °C
Analyte description	Neuron-specific enolase (NSE) has been detected in patients with certain tumors, namely: neuroblastoma, small cell lung cancer, medullary thyroid cancer, carcinoid tumors, pancreatic endocrine tumors, and melanoma. Studies of NSE as a tumor marker have concentrated primarily on patients with neuroblastoma and small cell lung cancer. Measurement of NSE levels in patients with these two diseases can provide information about the extent of the disease and the patient's prognosis, as well as about the patient's response to treatment.

PARAMETERS TESTED ON EACH LOT

Product appearance	Liquid, may turn slightly opaque during storage
Product concentration	5.0 mg/ml (+/- 10 %)
Immunoreactivity	80-120 % compared to the reference sample in an FIA test
IEF Profile	7.2 - 8.0
Purity	≥ 95 %

PARAMETERS DETERMINED DURING PRODUCT DEVELOPMENT

Subclass	IgG ₁
Association rate constant	Not Determined (N/D)
Dissociation rate constant	N/D
Affinity constant	2.2 x 10 ⁸ l/mol
Determination method	SPR analysis (Biacore)
Determination antigen	NSE, Scripps Laboratories (Cat N0224, Lot BF249001)

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB is a registered trademark of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.

Oy Medix Biochemica Ab

Headquarters and manufacturing
Klovinpellontie 3
FI-02180 Espoo, Finland

medix@medixbiochemica.com
www.medixbiochemica.com
VAT reg. no. FI4631532

Cross-reactivities

NNE (non-neuronal enolase) < 1.4 %

Epitope

Group B as described in Paus et al. (2011)

Pair recommendations

CAPTURE ANTIBODY

DETECTION ANTIBODY

9602
9601

9601
9602

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Product stability

TEMPERATURE, TIME

RESULT

-70 °C, 21 days
-20 °C, 21 days
+4 °C, 21 days
+35 °C, 7 days
+35 °C, 21 days
+45 °C, 3 days
+45 °C, 7 days

N/D
N/D
N/D
N/D
N/D
N/D
N/D

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2-8 °C in the product buffer.

Miscellaneous

Note that antibody binding requires Mg²⁺-ions and is sensitive to chelating agents (EDTA, EGTA, citrate).

References

Paus, E., Hirzel, K., Lidqvist, M., Höyhty, M., and Warren, D.J. (2011) TD-12 workshop report: characterization of monoclonal antibodies to neuron-specific enolase. Tumor Biol. 32:819-829.

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB is a registered trademark of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.

Oy Medix Biochemica Ab

Headquarters and manufacturing
Klovinpellontie 3
FI-02180 Espoo, Finland

medix@medixbiochemica.com
www.medixbiochemica.com
VAT reg. no. FI4631532