

Product Datasheet

Anti-Human CD21 Monoclonal Antibody, Azide Free Clone B-D6 854.170.005

Product Name	Anti-Human CD21 Monoclonal Antibody, Azide Free Clone B-D6
Catalog Number	854.170.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	ELISA
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Source	Daudi cell line
Specificity	Recognises the Epstein-Barr virus receptor (EBV R), CR2, C3d R, a 110 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells
Immunisation	Daudi cell line