

Product Datasheet

Plasma, Unit, Anti-TG Positive D1098

Product Name	Plasma, Unit, Anti-TG Positive
Catalog Number	D1098
Description	Each unit is approximately 800mL per donor though volumes can vary. Inquire for larger volumes, concentration ranges, testing, evaluation, or custom processing.
Tested Applications	QC Materials
Alternative Names	anti-TG, Plasma, Unit, Positive, thyroid, thyroglobulin, thyro globulin, TG, autoimmune, blood, auto immune, immune, immunology, blood, serology
Brand	Bioresource Technology
Form/Appearance	Liquid
Storage	-20°C
Source	Single Human Donor
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material, it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	Anticoagulant (typically sodium citrate)