

Product Datasheet

Mouse Immunoglobulin (IgG) Antibody JP204

Product Name	Mouse Immunoglobulin (IgG) Antibody
Catalog Number	JP204
Description	Affinity Purified Polyclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF, Blocking
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Polyclonal Antibody
Storage	-20 °C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 90 % by SDS-PAGE
Buffer	10mM Phosphate Buffer, 150mM Sodium Chloride, pH 7.4
Preservatives	0.1% ProClin 300