

Product Datasheet

Anti-h HBP 12202 PBS-5 101007

Product Name	Anti-h HBP 12202 PBS-5
Catalog Number	101007
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Species/Host	Mouse
Tested Applications	FIA
Alternative Names	Heparin binding protein, Azurocidin
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8 °C
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 95 %
Buffer	Phosphate-buffered saline (PBS) pH 7.4 with 0.095 % NaN ₃ as a preservative
IEF Profile	5.8-6.8
Specificity	Antibody recognizes human heparin binding protein (HBP)
Matched Pairs	Following pairs are especially recommended for the below mentioned assays: FIA: 12201 (capture) - 12202 (detection), and 12202 - 12201, 12201 - HM722, B-R37 - 12202
Shelf Life	12 months