

Product Datasheet

Anti-h IFN-g 13701 STPN-5 101008

Product Name	Anti-h IFN-g 13701 STPN-5
Catalog Number	101008
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Species/Host	Mouse
Tested Applications	LF, FIA, FC
Alternative Names	IFNG, IFN-?
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8 °C
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 95 %
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
IEF Profile	6.2-7.4
Specificity	Antibody recognizes human interferon gamma (IFN-g)
Matched Pairs	Following pairs are especially recommended for the below mentioned assays: FIA: 13701 (capture) - 13702 (detection), and 13702 - 13701
Shelf Life	Unspecified