

# Product Datasheet

## Urine, Single Donor, Mid-Stream Catch 991-03-S-MS

Product Name	Urine, Single Donor, Mid-Stream Catch
Catalog Number	991-03-S-MS
Description	Each sample is a minimum of 50mL per donor. Sample is collected during the mid-point of a urination event. Inquire for larger volumes per donor or custom criteria.
Tested Applications	QC Materials
Alternative Names	Urea, excretory fluid, pee, bodily waste, urinary, liquid waste, early, first, urinalysis
Brand	Lee Biosolutions
Form/Appearance	Liquid
Storage	-20°C
Source	Single Human Donor
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material, it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	None added, inquire if needed