

## Product Datasheet

### LowCross-Buffer® STRONG CR102

<b>Product Name</b>	LowCross-Buffer® STRONG
<b>Catalog Number</b>	CR102
<b>Description</b>	Antibody and sample diluent for minimizing nonspecific binding, cross-reactivities, and matrix effects in immunoassays
<b>Tested Applications</b>	LF, ELISA, WB, IHC, Immuno-PCR, Multiplex Assay
<b>Alternative Names</b>	interference reducing diluent
<b>Brand</b>	CANDOR
<b>Storage</b>	2 – 8 °C or -15 to -30 °C (tolerates repeated freezing and thawing cycles)
<b>Note</b>	LowCross-Buffer® is ready-to-use. Shake thoroughly before use
<b>Preservatives</b>	Contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1)