
Product Datasheet

LowCross-Buffer® STRONG CR102

Product Name	LowCross-Buffer® STRONG
Catalog Number	CR102
Description	Antibody and sample diluent for minimizing nonspecific binding, cross-reactivities, and matrix effects in immunoassays
Tested Applications	LF, ELISA, WB, IHC, Immuno-PCR, Multiplex Assay
Alternative Names	interference reducing diluent
Brand	CANDOR
Storage	2 - 8 °C or -15 to -30 °C (tolerates repeated freezing and thawing cycles)
Note	LowCross-Buffer® is ready-to-use. Shake thoroughly before use
Preservatives	Contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1)