
Product Datasheet

LowCross-Buffer® MILD protein-free CR301

| | |
|----------------------------|--|
| Product Name | LowCross-Buffer® MILD protein-free |
| Catalog Number | CR301 |
| Description | Protein- and animal-free antibody and sample diluent for minimizing nonspecific binding, cross-reactivities and matrix effects in immunoassays |
| Tested Applications | WB, IHC, Multiplex Assay, Immuno-PCR, ELISA, LF |
| Alternative Names | interference reducing diluent |
| Brand | CANDOR |
| Storage | 2 - 8 °C or -15 to -30 °C (tolerates repeated freezing and thawing cycles) |
| Preservatives | Contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1) |