
Product Datasheet

LowCross-Buffer® STRONG protein-free CR302

Product Name	LowCross-Buffer® STRONG protein-free
Catalog Number	CR302
Description	Protein- and animal-free antibody and sample diluent for minimizing nonspecific binding, cross-reactivities and matrix effects in immunoassays
Tested Applications	WB, IHC, Multiplex Assay, Immuno-PCR, ELISA, LF
Alternative Names	interference reducing diluent
Brand	CANDOR
Storage	2 - 8 °C or -15 to -30 °C (tolerates repeated freezing and thawing cycles)
Preservatives	Contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1)