

Product Datasheet

Sample Buffer CR105

Product Name	Sample Buffer
Catalog Number	CR105
Description	Sample and antibody dilution buffer for immunoassays
Tested Applications	WB, IHC, Multiplex Assay, Immuno-PCR, ELISA
Alternative Names	diluent
Brand	CANDOR
Storage	2 – 8 °C or -15 to -30 °C (tolerates repeated freezing and thawing cycles)
Note	Sample Buffer is ready-to-use. Shake thoroughly before use
Preservatives	Contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1)