

Product Datasheet

Sample Buffer for IHC CR106

Product Name	Sample Buffer for IHC
Catalog Number	CR106
Description	Sample and antibody dilution buffer for immunoassays
Tested Applications	IHC
Alternative Names	diluent for immunohistochemistry
Brand	CANDOR
Storage	2 – 8 °C or -15 to -30 °C (tolerates repeated freezing and thawing cycles)
Note	Sample Buffer is ready-to-use. Shake thoroughly before use
Preservatives	Contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1)