

Product Datasheet

The Blocking Solution CR110

Product Name	The Blocking Solution
Catalog Number	CR110
Description	Casein-based blocker for saturating free binding sites on plastic or other protein binding surfaces and for minimizing nonspecific binding to surfaces
Tested Applications	WB, IHC, Multiplex Assay, Immuno-PCR, ELISA, LF
Alternative Names	casein-based surface blocker
Brand	CANDOR
Storage	2 - 8 °C or -15 to -30 °C (tolerates repeated freezing and thawing cycles)
Note	The Blocking Solution is ready-to-use. Shake thoroughly before use
Preservatives	Contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1)