
Product Datasheet

SmartBlock™ CR113

Product Name	SmartBlock™
Catalog Number	CR113
Description	BSA-free blocker for saturating free binding sites on plastic surfaces or other protein binding surfaces
Tested Applications	WB, Multiplex Assay, Immuno-PCR, ELISA, LF
Alternative Names	peptide-based surface blocker
Brand	CANDOR
Storage	2 - 8 °C or -15 to -30 °C (tolerates repeated freezing and thawing cycles)
Note	SmartBlock™ is ready-to-use. Shake thoroughly before use
Preservatives	Contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1)