

Product Datasheet

Liquid Plate Sealer® animal-free CR163

Product Name	Liquid Plate Sealer® animal-free
Catalog Number	CR163
Description	Stabilizer for coated antibodies and antigens on polystyrene- or glass surfaces. Free from animal derived ingredients and animal proteins
Tested Applications	LF, ELISA, Immuno-PCR, Multiplex Assay
Alternative Names	coating stabilizer
Brand	CANDOR
Storage	2 – 8 °C
Note	Liquid Plate Sealer® animal-free is ready-to-use. Liquid Plate Sealer® animal-free can either be added to the finished coated microtiter plates, polystyrene beads, glass slides etc. (step A) or alternatively added directly to
Preservatives	Contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1)