

## Product Datasheet

### Liquid Plate Sealer® animal-free CR163

<b>Product Name</b>	Liquid Plate Sealer® animal-free
<b>Catalog Number</b>	CR163
<b>Description</b>	Stabilizer for coated antibodies and antigens on polystyrene- or glass surfaces. Free from animal derived ingredients and animal proteins
<b>Tested Applications</b>	Multiplex Assay, Immuno-PCR, LF, ELISA
<b>Alternative Names</b>	coating stabilizer
<b>Brand</b>	CANDOR
<b>Storage</b>	2 – 8 °C
<b>Note</b>	Liquid Plate Sealer® animal-free is ready-to-use. Liquid Plate Sealer® animal-free can either be added to the finished coated microtiter plates, polystyrene beads, glass slides etc. (step A) or alternatively added directly to
<b>Preservatives</b>	Contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1)