

Product Datasheet

Assay Defender® CR180

Product Name	Assay Defender®
Catalog Number	CR180
Description	Sample diluent for blocking HAMA and other high affinity interfering antibodies and minimizing nonspecific binding, cross-reactivities and matrix effects in immunoassays based on human or animal body fluids
Tested Applications	Multiplex Assay, Immuno-PCR, ELISA
Alternative Names	HAMA and interference blocker
Brand	CANDOR
Storage	2 – 8 °C
Preservatives	Contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1)