

Product Datasheet

Sample Buffer protein-free CR305

Product Name	Sample Buffer protein-free
Catalog Number	CR305
Description	Protein- and animal-free sample and antibody dilution buffer for immunoassays
Tested Applications	WB, IHC, Multiplex Assay, Immuno-PCR, ELISA
Alternative Names	diluent
Brand	CANDOR
Storage	2 – 8 °C or -15 to -30 °C (tolerates repeated freezing and thawing cycles)
Preservatives	Contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1)