

Product Datasheet

HRP-Protector™ protein-free CR322

Product Name	HRP-Protector™ protein-free
Catalog Number	CR322
Description	Protein- and animal-free diluent for long-term storage of HRP conjugates
Tested Applications	WB, IHC, Multiplex Assay, ELISA
Alternative Names	HRP-conjugate Stabilizer
Brand	CANDOR
Storage	2 – 8 °C (Does not tolerate freezing)
Preservatives	Contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1)