

Product Datasheet

Antibody Stabilizer TRIS protein-free CR330

Product Name	Antibody Stabilizer TRIS protein-free
Catalog Number	CR330
Description	Protein- and animal-free stabilizer for long-term storage of proteins or antibodies at 2 – 8 °C
Tested Applications	WB, IHC, Multiplex Assay, Immuno-PCR, ELISA
Alternative Names	stabilizer
Brand	CANDOR
Storage	2 – 8 °C (Does not tolerate freezing)
Preservatives	Contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1)