

Product Datasheet

FSH Antigen, Recombinant 610174

Product Name	FSH Antigen, Recombinant
Catalog Number	610174
Description	Recombinant Human Follicle Stimulating Hormone (FSH) protein, including both alpha and beta subunits, with a C-terminal histidine tag in the beta subunit. Predicted molecular weight: 24 kDa
Species/Host	CHO Cells
Alternative Names	FSH, Follicle-Stimulating Hormone
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Storage	2-8°C
Purity	≥ 90% (HPLC-SEC)
Source	CHO cells
Associated Products	Anti-h FSH 6602 and FA311, Anti-h Alpha subunit 5501, 5503 and 6601
Buffer	Lyophilized from 50 mM Tris-HCl, 150 mM NaCl, pH 7.5
Molecular Weight	~ 24 kDa
Shelf Life	Unspecified
Biological Activity	≥ 3700 IU/mg (immunoassay)
Unit Definition	1 IU is defined as the amount of FSH that has an activity corresponding to 113.88 ug of pure human urinary FSH
Tag	His-tag