

Product Datasheet

Cotinine-4 Antibody PP301IF

Product Name	Cotinine-4 Antibody
Catalog Number	PP301IF
Species/Host	Rabbit
Tested Applications	ELISA, LF
Immunogen	Cotinine-4-BgG
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Polyclonal Antibody, Ig Fraction
Storage	Short Term (\leq 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	Phosphate Buffered Saline, pH 7.0
Preservatives	0.05% Sodium Azide