

# Product Datasheet

## Testosterone Depleted Plasma H5102

|                     |                                                                        |
|---------------------|------------------------------------------------------------------------|
| Product Name        | Testosterone Depleted Plasma                                           |
| Catalog Number      | H5102                                                                  |
| Description         | Testosterone Depleted, Pooled, Processed, & 0.2 Micron Filtered Plasma |
| Tested Applications | QC Materials, Base Matrix                                              |
| Brand               | Bioresource Technology                                                 |
| Form/Appearance     | Liquid                                                                 |
| Storage             | -20°C                                                                  |
| Source              | Pooled Human Donors                                                    |