

Product Datasheet

Anti-Human CD28 Monoclonal Antibody, PE Conjugated Clone B-T3 854.222.010

Product Name	Anti-Human CD28 Monoclonal Antibody, PE Conjugated Clone B-T3
Catalog Number	854.222.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	PE
Tested Applications	FC
Immunogen	U266 cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1% BSA and 0.1% sodium azide
Storage	Stable at +2-8°C for six months after reconstitution. DO NOT FREEZE.
Note	Reconstitute with 1 ml deionised water Use 10 µl to label 106 cells or 100 µl of whole blood
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the B7/BB1 receptor (B7/BB1 R), a 44 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells