

Product Datasheet

Pregnancy Associated Plasma Protein A, (PAPP-A) Placental Blood 496-91

Product Name	Pregnancy Associated Plasma Protein A, (PAPP-A) Placental Blood
Catalog Number	496-91
Tested Applications	Control, Calibration
Brand	Lee Biosolutions
Form/Appearance	Clear to opaque, colorless to light yellow frozen liquid
Concentration	0.5 to 1.5 mg/mL protein (Coomassie)
Storage	-20°C
Purity	≥ 90% (SDS-PAGE)
Source	Human Placental Blood
Buffer	Phosphate buffered saline, pH 7.4
Shelf Life	1 year