

Product Datasheet

Chikungunya virus Antigen, recombinant LA830

Product Name	Chikungunya virus Antigen, recombinant
Catalog Number	LA830
Brand	EastCoast Bio
Form/Appearance	Liquid; Purified Recombinant Antigen
Storage	-20°C. Avoid repeated freezing and thawing.
Epitope	E1 protein
Purity	≥95%
Buffer	DPBS, pH 7.4 with 5 mM EDTA, 1 µg/ml Leupeptin, Aprotinin and Pepstatin A
Molecular Weight	50 kDa
Expression System	Insect Cell
Preservatives	0.1% Thimerosal