

Product Datasheet

Chikungunya virus Antigen, recombinant LA825

Product Name	Chikungunya virus Antigen, recombinant
Catalog Number	LA825
Species/Host	E. coli
Brand	EastCoast Bio
Form/Appearance	Liquid; Purified Recombinant Antigen
Storage	-20°C. Avoid repeated freezing and thawing.
Epitope	Capsid protein
Purity	≥90%
Buffer	20 mM Phosphate Buffer, pH 7.4 with 0.01% SDS
Molecular Weight	35.8 kDa