

Product Datasheet

Anti-HIV-1 p24 R14102 PBS-5 140072

| | |
|----------------------------|--|
| Product Name | Anti-HIV-1 p24 R14102 PBS-5 |
| Catalog Number | 140072 |
| Description | Recombinant rabbit antibody, cultured in vitro under conditions free from animal-derived components. |
| Species/Host | Rabbit |
| Tested Applications | WB, ELISA, LF, FIA |
| Immunogen | Human immunodeficiency virus 1 |
| Brand | Medix Biochemica |
| Form/Appearance | Liquid, may turn slightly opaque during storage |
| Concentration | 5.0 mg/ml (+/-10 %) |
| Storage | +2-8 °C |
| Isotype | Rabbit IgG |
| Clonality | Monoclonal |
| Epitope | Not determined |
| Purity | ≥ 95 % |
| Buffer | Phosphate-buffered saline (PBS) pH 7.4 with 0.095 % NaN3 as a preservative |
| Cross Reactivity | Recognizes HIV-1 p24 group O and N. Antibody also recognizes HIV-1 p24 group M clades A, B, C and D. |
| Specificity | Antibody recognizes human immunodeficiency virus 1 p24 (HIV-1 p24) |
| Matched Pairs | FIA: R14102 (capture) - R14101 (detection) LF: R14102 (membrane) - R14101 (particles) |
| Shelf Life | Unspecified |