

Product Datasheet

Treponema pallidum TpN 15,17,47 antigen, recombinant
610182

Product Name	Treponema pallidum TpN 15,17,47 antigen, recombinant
Catalog Number	610182
Description	Recombinant Treponema pallidum protein includes sequences of TpN 15, 17 and 47 with a C-terminal histidine tag
Brand	Medix Biochemica
Form/Appearance	Liquid
Concentration	≥ 1.0 mg/ml
Storage	−20 °C
Purity	≥ 80% (CE-SDS)
Buffer	50 mM Tris-HCl, 150 mM NaCl, 1.0 mM Tris(2-carboxyethyl)phosphine, 0.02% (v/v) Tween 20
Molecular Weight	Predicted molecular weight: 75 kDa
Matched Pairs	V2101
Shelf Life	Unspecified
Expression System	E. coli