

Product Datasheet

Rheumatoid Factor Concentrate; ≥550 IU/mL
508-01

Product Name	Rheumatoid Factor Concentrate; ≥550 IU/mL
Catalog Number	508-01
Tested Applications	Blocking
Brand	Lee Biosolutions
Form/Appearance	Liquid in defibrinated plasma
Concentration	Reported in mg/mL (Biuret)
Storage	-20°C
Source	Defibrinated human plasma
Shelf Life	1 year from date of manufacture
Biological Activity	≥ 550 IU/mL (Turbidimetric Immunoassay)
Specific Activity	Reported in IU/mg protein