

Product Datasheet

Anti-Human CD34 Class I Monoclonal Antibody, Azide Free Clone B-G25 854.240.005

Product Name	Anti-Human CD34 Class I Monoclonal Antibody, Azide Free Clone B-G25
Catalog Number	854.240.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	Functional assay, FC
Immunogen	KG-1a cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Gp105-120 antigen, a 105-120 kDa protein
Biological Activity	Triggers KG-1a cell aggregation
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells