

Product Datasheet

LowCross-Buffer(R) (5x) CR303

Product Name	LowCross-Buffer(R) (5x)
Catalog Number	CR303
Description	Protein- and animal-free ready-to-use solution for the storage of AP-conjugates at 2 - 8 °C. To be diluted into an exisiting sample dilution buffer or into distilled water.
Tested Applications	WB, IHC, Multiplex Assay, Immuno-PCR, LF, ELISA
Brand	CANDOR
Form/Appeara	Clear liquid
Storage	2 - 8 °C
Associated Products	CR100, CR300