

## Product Datasheet

### LowCross-Buffer(R) (5x) CR303

<b>Product Name</b>	LowCross-Buffer(R) (5x)
<b>Catalog Number</b>	CR303
<b>Description</b>	Protein- and animal-free ready-to-use solution for the storage of AP-conjugates at 2 - 8 °C. To be diluted into an existing sample dilution buffer or into distilled water.
<b>Tested Applications</b>	WB, IHC, Multiplex Assay, Immuno-PCR, ELISA, LF
<b>Brand</b>	CANDOR
<b>Form/Appearance</b>	Clear liquid
<b>Storage</b>	2 - 8 °C
<b>Associated Products</b>	CR100, CR300