

Product Datasheet

Nitazene-BSA Antigen LA824

Product Name	Nitazene-BSA Antigen
Catalog Number	LA824
Description	Nitazene-BSA Antigen
Tested Applications	ELISA, LF
Brand	Medix Biochemica
Form/Appearance	Liquid, Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	>90%
Buffer	10 mM Phosphate buffer, pH 7.2, 0.01% sodium azide as a preservative
Matched Pairs	HM1453 (detection)
Shelf Life	Unspecified