

Product Datasheet

Nitazene Antibody HM1453

Product Name	Nitazene Antibody
Catalog Number	HM1453
Description	Mouse anti-Isotonitazene /Nitazene Monoclonal Antibody
Tested Applications	ELISA, LF
Immunogen	Isotonitazene-KLH
Brand	Medix Biochemica
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	Isotonitazene 100%, Butonitazene 80%, Etonitazene 100%, Etodesnitazene 80%, Metonitazene 90%, Metodesnitazene 100%, N-Pyrrolidino Etonitazene 100% No cross reactivity with Amphetamine, Methamphetamine, THC, Opiates, TCA, PCP, Methadone, MDMA, Oxycodone or LSD
Specificity	Nitazene
Matched Pairs	Nitazene-BSA Antigen LA824 (capture)
Shelf Life	Unspecified